

Capitainer® B50

Quantitative blood sampling card

Next generation microsampling card

Capitainer®B50 is the blood sampling card for a larger sample volume, enabling microsampling for a wider range of analytes and application techniques. All Capitainer®B products deliver volumetric accuracy and precision with our technology.

It empowers healthcare providers, researchers and laboratories across the world to provide accurate and easy to use at-home sampling solutions for their customers and patients.



Accurate

Automatically meters
2×50µl of blood



Reliable

Designed to minimize
the risk of over- or
under-fill



Trusted

Successful sampling
indicator integrated in
the device



Protected sample

Sample protected from
contamination within the device



User-friendly

Allows applications of
multiple drops until
metering channel is filled