

Capitainer® DIP70

Volumetric urine sampling card

Next generation microsampling card

Capitainer®DIP70 is developed to facilitate remote urine sampling. The product collects and stores two volume defined samples in dried format. The drying of the sample facilitates shipping via regular mail reducing the risk of microbial growth and degradation of analytes.

It empowers healthcare providers, researchers and laboratories across the world to provide accurate and easy to use at-home urine sampling solutions for their customers and patients.



Accurate

Automatically meters
2×70µl of urine



User-friendly

Easy to use — simply dip
to collect the sample



Hygienic

No need to handle
urine tubes



Shipping

Fits standard envelopes for
convenient and cost-effective
transport



Stable

Dried format minimizes
sample contamination
by unwanted fungal
growth during postal
transport