

Capitainer®SEP10

Quantitative blood separation card

The Revolution has Started

Capitainer®SEP10 is developed for collection of blood where the cells are separated from the blood at the site of collection without centrifugation. The cell separated blood matrix is metered to an exact volume and dried for best shipping stability in ambient temperature.

Capitainer®SEP10 generates a plasma-like sample ideal for testing of substances that cannot be measured from a whole blood sample.



Accurate

Automatically generates
2×10 μ L cell separated blood



Reliable

Microfluidic control
of correct volume



Trusted

Successful indicator
integrated in the device



Protected sample

Sample protected from
contamination



User-friendly

Separated blood without
centrifugation